

## **CO<sub>2</sub> exchange in plants**

### **The effect of light intensity and quality on photosynthetic CO<sub>2</sub> fixation and respiration in C3 and C4 plants**

Using light energy and chlorophyll plants are able to oxidize water and use these electrons to generate energy and to reduce carbon that is taken up from the atmosphere in the form of CO<sub>2</sub>. This is the process of photosynthesis and the end products are energy-rich sugar molecules. The overall stoichiometry of photosynthesis has one mole (or molecule) of oxygen is produced per mole of CO<sub>2</sub> that is fixed or assimilated. This allows us to measure the apparent photosynthetic rate by assessing the amount of O<sub>2</sub> that a leaf produces or the amount of CO<sub>2</sub> that is taken up by a leaf. However, these measurements are complicated by the fact that plants are also respiring as they release the energy stored in sugars and other energy-rich molecules; this process consumes O<sub>2</sub> and releases CO<sub>2</sub> thereby effectively working opposite to what is achieved during photosynthesis.

When we measure net CO<sub>2</sub> uptake by a leaf we do not distinguish between photosynthetic and respiratory gas exchange and simply measure the sum of these two processes. To know the true rate of photosynthesis we need to be able to measure leaf CO<sub>2</sub> output in the light so as to know the rate of respiration. However, simultaneously measuring photosynthetic CO<sub>2</sub> uptake and respiratory CO<sub>2</sub> output is not easy and requires the use of <sup>14</sup>C-labeled substrates or tritiated H<sub>2</sub>O or <sup>18</sup>O<sub>2</sub>. Furthermore, the rate of CO<sub>2</sub> uptake at saturating light intensities may be 30 times the respiratory rate but this situation is complicated by the impact of photorespiration (light-dependent CO<sub>2</sub> output), which varies from plant to plant and is strongly influenced by whether the plant uses C4 or C3 photosynthesis. In C3 plants rates of photorespiration can account for a sizeable fraction of the total gas exchange whereas in C4 plants photorespiration rates are negligible.

A convenient and widely used method to measure the apparent photosynthesis rate (APS) is to use an infrared gas analyzer (IRGA). CO<sub>2</sub> absorbs strongly in certain infrared wavelengths and this absorption is used to measure the CO<sub>2</sub> concentration in a volume of gas. The principle behind the IRGA is to measure the IR absorption and to convert this absorption into an electrical signal and thence to a concentration of CO<sub>2</sub> in parts per million (ppm).

### **Objectives:**

During this exercise you will learn:

- How to use an IRGA to measure changes in CO<sub>2</sub> concentrations
- The effects of light intensity on the rate of CO<sub>2</sub> exchange in two different type of plants (C3 and C4) and how to determine the light compensation point from a light response curve.
- The effects of O<sub>2</sub> concentration on the rate of CO<sub>2</sub> exchange in C3 and C4 plants

- How to use a radiometer and the quantum meter to measure the light intensities and photosynthetically active radiation (PAR) or photosynthetic photon flux density (PPFD), respectively.

## **Materials**

You are provided with:

Plant material (bean and corn)  
Calibrated IRGA and light source  
Gas cylinder filled with air and/or other gas mixes  
Acetate grid to measure leaf area  
LabPro connected to a computer

### **Experiment one: using the IRGA to measure apparent photosynthesis**

Procedure:

- 1) Set up the IRGA apparatus as follows: Connect the pump to the flow meter to the inlet of the leaf chamber. Connect the outlet of the leaf chamber to the humidity sensor to the drying column and finally to the IRGA.
- 2) Attach a gas bag containing air to the inlet of the pump. Let the system flush for 3 – 5 minutes and record the  $[\text{CO}_2]$  as reference  $[\text{CO}_2]$ .
- 3) Set the Logger Pro display so that the y-axis of the graph has a range of approx. 130 ppm  $\text{CO}_2$ . This range should include values that are up to 30 ppm above the reference  $\text{CO}_2$  concentration and up to 100 ppm below the reference value.
- 4) Seal a leaf inside the leaf chamber and place the LED light source on top of the leaf chamber
- 5) The leaf chamber has four gas ports that are grouped in pairs. Each pair has a port on the upper surface of the chamber that is located directly above a port on the lower surface of the chamber. These ports distribute the gases to the upper and lower leaf surfaces through manifolds. Each pair of ports is attached to plastic tubing that connects via a Y-connector to a single piece of tubing.
- 6) Turn on the light adjusting to the maximum output.
- 7) You will observe a decline in the  $\text{CO}_2$  concentration of the gas leaving the leaf chamber as photosynthesis consumes  $\text{CO}_2$  delivered in the reference gas. This may take several minutes.
- 8) When a steady value of  $\text{CO}_2$  is reached record the  $\text{CO}_2$  value as your “Analysis  $\text{CO}_2$ ” value.

- 9) Turn the light source off. The CO<sub>2</sub> concentration will increase as photosynthesis stops. You should see a peak in the CO<sub>2</sub> level within the first few minutes. Record the CO<sub>2</sub> concentration attained after the concentration has reached a steady level. This is your “dark analyses CO<sub>2</sub>” value.
- 10) Click the STOP button on the computer screen to stop data collection and save your data (“File”, “Save As”) giving it a name.
- 11) Measure the area of your leaf. Note: if the leaf completely filled the chamber the enclosed area is 9 cm<sup>2</sup>. If your leaf did not completely fill the chamber you must use the acetate grid to estimate leaf area. Place the grid on the surface of the chamber.
- 12) Using the acetate grid count the number of interstices that are completely enclosed by the area of the leaf and assign a value of 1 to each. Count the number of interstices that fall on the leaf margin and assign these a value of 0.5. Total the number of interstices and divide that value by 4 to obtain the leaf area in cm<sup>2</sup>.
- 13) Remove the leaf from the leaf chamber. KEEP your leaf we may need it to determine chlorophyll content/ fresh weight/ dry weight if there is enough time.
- 14) Repeat for the corn plant. In addition, there are other plants available that have been exposed to a variety of environmental stresses. Measure the CO<sub>2</sub> uptake rate for one of these. The plant you have just measured is a control, unstressed plant.

### Calculation of the CO<sub>2</sub> exchange rate

Measurements of CO<sub>2</sub> exchange rates in leaves are typically expressed in units of CO<sub>2</sub> exchange per unit time per unit leaf area (μmol CO<sub>2</sub>/m<sup>2</sup>/s). There are other ways that we can express these data such as CO<sub>2</sub> exchange based on leaf weight (dry or fresh) or leaf chlorophyll content. To express your data as μmol CO<sub>2</sub>/m<sup>2</sup>/s you need to calculate the following:

- 1) The difference between the CO<sub>2</sub> concentration in the reference and analyses gases. i.e.  $\Delta\text{CO}_2 = \text{Reference CO}_2 - \text{analysis CO}_2$
- 2) Convert the  $\Delta\text{CO}_2$  value from (ppm) to (μmol CO<sub>2</sub>/L) using the following calculation:

$$\Delta\text{CO}_2 / 22.415 \times ([T+C]/T)$$

C = temperature in °C

T = absolute temperature (273 K)

- 3) Multiply the  $\Delta\text{CO}_2$  by the flow rate (L/s) to get the CO<sub>2</sub> exchange rate per second.
- 4) Divide the CO<sub>2</sub> exchange rate by the leaf area (m<sup>2</sup>) to get μmol CO<sub>2</sub>/m<sup>2</sup>/s.

## Experiment 2. Measuring the light dependence of photosynthesis

Procedure:

- 1) Set up the IRGA and proceed as for experiment one.
- 2) Seal a leaf from a bean plant inside the leaf chamber.
- 3) Turn on the light source and adjust to 10% of the maximum output. Record the light output.
- 4) Observe the decline in CO<sub>2</sub> concentration in the gas leaving the leaf chamber. Record the CO<sub>2</sub> concentration at steady state conditions. This is the “analysis CO<sub>2</sub>” at 10% light output.
- 5) Increase the light output to 20% of the maximum level. Measure the light level.
- 6) Observe the change in the CO<sub>2</sub> concentration of the gas leaving the leaf chamber. When the CO<sub>2</sub> concentration reaches a new steady state level record the concentration.
- 7) Increase the light output to 40% and repeat as for steps 5) and 6). Repeat this for a light output of 60% and 80% of the maximum output and for the maximum light level. Remember to measure the light level in each case.
- 8) Turn off the light and observe the increase in CO<sub>2</sub> concentration. Record the CO<sub>2</sub> concentration when it reaches a steady state level.
- 9) Stop data collection and follow steps 10) to 13) as for experiment one. Calculate the CO<sub>2</sub> exchange rates.
- 10) IF there is enough time you can repeat this experiment for the corn plant.
- 11) Record your results in the table below

[illegible]

Construct a photosynthetic light response curve for the bean and the corn plant.

Determine:

- The light compensation point
- The rate of dark respiration
- The photochemical efficiency
- The light saturation point

See the end of the handout for a light response curve and information on these parameters.

### **Experiment 3. Photosynthetic response to O<sub>2</sub> concentration in a C3 versus a C4 plant**

Procedure:

- 1) Set up the apparatus as in step 1) of experiment one with the exception that you will attach a high O<sub>2</sub> air bag to the flow meter.
- 2) Allow the system to flush for 3 minutes.
- 3) **Record the CO<sub>2</sub> concentration as your “reference CO<sub>2</sub> concentration- high O<sub>2</sub>”.**
- 4) Set the Logger Pro display so that the y-axis of the graph has a range of approx. 130 ppm CO<sub>2</sub>. This range should include values that are up to 30 ppm above the reference CO<sub>2</sub> concentration and up to 100 ppm below the reference value.
- 5) Seal a bean leaf inside the leaf chamber as for step 2) of experiment two and proceed to step 8).
- 6) Stop data collection by clicking on the STOP button on the computer screen and save your data.
- 7) Detach the high O<sub>2</sub> air bag from the inlet of the pump and seal it. Do not remove the leaf!
- 8) Attach a low O<sub>2</sub> air bag to the flow meter. Flush the system for 3 minutes and then record the CO<sub>2</sub> concentration as your “reference CO<sub>2</sub> concentration – low O<sub>2</sub>”. Repeat steps 4) to 6). You will omit step 5) because you have not removed your leaf.
- 9) Measure the leaf area using the acetate grid and remove the leaf.
- 10) Repeat for corn.
- 11) Record your results in the table:



**Light Compensation point** is the light level at which CO<sub>2</sub> consumption by photosynthesis exactly matches CO<sub>2</sub> production arising from respiration. Extrapolate the linear portion of your light response curve to intercept the *x-axis* at the point where the photosynthesis rate is zero.

**Rate of dark respiration** can be estimated from your light response curve by extrapolating it to intercept the *y-axis* at the point where the light intensity is zero.

**Photochemical efficiency** is obtained by calculating the initial slope of the light response curve. This value reflects the efficiency by which light energy is used for photosynthesis.

**Light saturation point of photosynthesis** is the light intensity at which the light response curve plateaus. Beyond this point light does not limit photosynthesis but other factors become limiting.

**A light response curve** (taken from “Plant Physiology” Taiz and Zeiger, 3<sup>rd</sup> edition).

